



EasySC Exosome DNA Purification Kit

Catalog Number: D035-1, D035-2

Table 1. Kit Components and Storage

Kit Component	D035-1 (25 preps)	D035-2 (50 preps)	Storage	Stability
Buffer AL	8 mL	15 mL	RT	The product is stable for one year when stored as directed.
Buffer DW1*	6 mL	15 mL	RT	
Buffer DW2*	10 mL	15 mL	RT	
Buffer AE	10 mL	15 mL	RT	
Mini Column	25	50	RT	
2 mL Collection Tube	25	50	RT	

* Prior to use, add absolute ethanol to **Buffer DW1** and **Buffer DW2** according to the bottle label.

Product Description

EasySC Exosome DNA Purification Kit is designed for the isolation of exosomal DNA from pure exosome isolated by our Exosome Isolation Kits (Cat. #: D030, D031, D032, D033, D034). The kit also features a buffer system that facilitates complete exosome lysis for efficient nucleic acid isolation. Exosomal DNA is bound to spin column, washed and eluted. Normally, 100-300 ng of exosomal DNA can be purified from the exosome isolated from 200 μ L blood serum or 5 mL cell media. The extracted DNA is ready for downstream applications such as sequencing, and PCR. No need to precipitate, concentrate or desalt.

Features

- ❖ Fast: The purification process takes only 30 minutes.
- ❖ Fully compatible: The isolated DNA can be used in most downstream applications such as PCR and sequencing.
- ❖ Clean: Minimal carryover of co-precipitating proteins.

Purification Protocol

1. Exosome isolation

Cell culture media: using our ExoFast™ Exosome Isolation Reagent from cell culture media (Cat. No. D030) to isolate exosome.

Serum: using our ExoFast™ Exosome Isolation Reagent from serum (Cat. No. D031) to isolate exosome.

Plasma: using our ExoFast™ Exosome Isolation Reagent from plasma (Cat. No. D032) to isolate exosome.

Urine: using our ExoFast™ Exosome Isolation Reagent from urine (Cat. No. D033) to isolate exosome.

Other body fluids (cerebrospinal fluid (CSF), ascitic fluid, amniotic fluid, milk, and saliva): using our ExoFast™ Exosome Isolation Reagent from other body fluids (Cat. No. D034) to isolate exosome.

2. Transfer 250 μ L of exosome suspension to a 1.5 mL centrifuge tube. If the volume of the sample is less than 250 μ L, bring the volume up to 250 μ L with PBS or TE buffer.
3. Add 250 μ L Buffer AL to the sample. Vortex to mix well, and incubate at 65°C for 15-30 min.
Note: Buffer AL may be precipitate during storage, if happen, heat it at 55°C to dissolve. If RNA need be removed, add 10 μ L RNase A Solution (10 mg/mL) to the sample.
4. Add 250 μ L isopropanol to the sample, and vortex for 30 second. Centrifuge briefly to collect any drops from the inside of the lid.
5. Insert a DNA Mini Column into a 2 mL Collection Tube. Transfer the samples from step 4 to the DNA Mini Column, then centrifuge at 10,000 x g for 1 min.
6. Discard the filtrate and reuse the collection tube. Add 600 μ L of Buffer DW1 to the DNA Mini Column, invert and mix once, then centrifuge at 10,000 x g for 1 min.
Note: Buffer DW1 must be diluted with absolute ethanol according to the bottle label before use.
7. Discard the filtrate and reuse the collection tube. Add 600 μ L of Buffer DW2 to the DNA Mini Column, then centrifuge at 10,000 x g for 1 min.
Note: Buffer DW2 must be diluted with absolute ethanol according to the bottle label before use.
8. Discard the filtrate and reuse the collection tube. Add 600 μ L of Buffer DW2 to the DNA Mini Column, then centrifuge at 10,000 x g for 1 min.
9. Discard the filtrate and reuse the collection tube. Centrifuge the empty DNA Mini Column at 12,000 x g for 3 min.
Note: This step is critical for removing of trace ethanol that may interfere with downstream applications.
10. Transfer the DNA Mini Column into a new nuclease-free 1.5 mL microcentrifuge tube, add 30-100 μ L Buffer AE preheated to 70°C. Let sit at RT for 3 min, then centrifuge at 10,000 x g for 1 min.
Note: To improve the yield, repeat this step for a second elution step.
11. Discard the column and store the DNA at -20°C.

Troubleshooting

Problem	Possible cause and suggestions
Column is clogged	<p>Incomplete sample lysis</p> <ul style="list-style-type: none"> • The sample was not mixed well with Buffer AL. After adding Buffer AL, invert and mix for 3-5 times, then vortex at maximum speed to mix the sample with Buffer AL. • Add too much sample. Reduce sample amount.
Low purity	<p>Incomplete sample lysis</p> <ul style="list-style-type: none"> • The sample was not mixed well with Buffer AL. After adding Buffer AL, invert and mix for 3-5 times, then vortex at maximum speed to mix the sample with Buffer AL. <p>Insufficient washing procedure</p> <ul style="list-style-type: none"> • The wash buffer Buffer DW1 and Buffer DW2 must be diluted with absolute ethanol before use.
Poor yield	<p>Incomplete sample lysis</p> <ul style="list-style-type: none"> • The sample was not mixed well with Buffer AL. After adding Buffer AL, invert and mix for 3-5 times, then vortex at maximum speed to



7040 Virginia Manor Road
Beltsville, MD 20705, USA
Web: www.abpbio.com
Email: info@abpbio.com

	<p>mix the sample with Buffer AL.</p> <p>Insufficient elution buffer volume</p> <ul style="list-style-type: none">• Increase elution buffer volume, and repeat elution step. <p>Wash buffer no ethanol added</p> <ul style="list-style-type: none">• The wash buffer Buffer DW1 and Buffer DW2 must be diluted with absolute ethanol before use.
--	--