



## CellCycle™ Green Stain

Catalog Number: A054

Table 1. Product Package and Storage

Material	Amount	Concentration	Storage	Stability
CellCycle™ Green Stain	100 µL	500X	-20 °C	One year when stored as directed.

Number of assays: 100 assays.

Approximate fluorescence excitation/emission maxima, in nm: 488/530, bound to DNA.

### Introduction

Live cell studies of cellular DNA content and cell cycle distribution are useful to detect variations of growth patterns due to a variety of physical, chemical, or biological means, to monitor apoptosis, and to study tumor behavior and suppressor gene mechanisms. In a given population, cells are distributed among three major phases of cell cycle: G<sub>0</sub>/G<sub>1</sub> phase (one set of paired chromosomes per cell), S phase (DNA synthesis with variable amount of DNA), and G<sub>2</sub>/M phase (two sets of paired chromosomes per cell, prior to cell division). DNA content can be measured using fluorescent, DNA-selective stains that exhibit emission signals proportional to DNA mass. Flow cytometric analysis of these stained populations is then used to produce a frequency histogram that reveals the various cell cycle phases. This analysis is typically performed on permeabilized or fixed cells using a cell-impermeant nucleic acid stain, but is also possible using live cells and a cell-permeant nucleic acid stain.

The CellCycle™ Green stain is DNA-selective, cell membrane-permeant, and nonfluorescent stain for DNA content analysis in living cells. The CellCycle™ Green stain is fluorescent upon binding to double-stranded DNA. CellCycle™ Green stain is excited using 488 nm laser with emission ~530 nm.

The staining protocol is simple and includes incubating suspended cells in the presence of CellCycle™ stain and directly measuring the fluorescence without the need for any additional treatment or centrifugation steps. This live cell stain allows the simultaneous co-staining of the cell population for other parameters, and allows for the possibility of cell sorting based on DNA content.

### Materials required but not provided

- Cells and culture medium
- Flow cytometer tubes

### Experimental Protocols

**Note:** The following staining protocol was optimized using Jurkat cells, a human T-cell leukemia line, in complete RPMI medium containing 10% fetal bovine serum with staining at 37°C, but can be adapted to most cell types. These stains can also be used for cells suspended in Hanks' Balanced Salt Solution (HBSS) or phosphate-buffered saline (PBS).

1. Remove the CellCycle™ Green stain from the freezer and allow it to equilibrate to room temperature.
2. Prepare flow cytometry tubes each containing 0.5 mL of cell suspension in complete media at a concentration of  $1 \times 10^6$  cells/mL.
3. To each tube add 1 µL of CellCycle™ Green stain.
4. Incubate at 37°C for 30 minutes, protected from light.
5. Analyze the samples without washing, using 488 nm excitation, and collect emission using a 530/30 bandpass filter or equivalent.