



7040 Virginia Manor Road
Beltsville, MD 20705, USA
Web: www.abpbio.com;
Email: info@abpbio.com

Product Information

Rhodamine 123

Cat. No.	Product Name	Unit Size
C044	Rhodamine 123	25 mg

Storage upon receipt:

- -20°C
- Protect from light

Product Description

Rhodamine 123 is a cell-permeant, cationic, green-fluorescent dye that is readily sequestered by active mitochondria without cytotoxic effects. This product has been used to assay mitochondrial membrane potential in populations of apoptotic cells..

Specifications:

Chemical name: Rhodamine 123

Molecular Formula: $C_{21}H_{17}ClN_2O_3$

Molecular Weight: 380.83

CAS Number: 62669-70-9