



Product Information

CDCFDA, SE

Cat. No.	Product Name	Unit Size
C036	CDCFDA, SE	25 mg

Storage upon receipt:

- -20°C
- Protect from light

Product Description

5(6)-CDCFDA, SE is a fixable-cell-permeant, fluorescein-based tracer for very long-term cell labeling.

Specifications:

Chemical name: 5-(and-6)-Carboxy-2',7'-dichlorofluorescein Diacetate, succinimidyl ester

Molecular Formula: C₂₅H₁₇Cl₂NO₁₁

Molecular Weight: 626.36

CAS Number: 147265-60-9