



405 E Gude Dr, STE 214
Rockville, MD 20850, USA
www.abpbio.com

CERTIFICATE OF ANALYSIS

Catalog Number C001-2
Product Name DAPI, 5 mg/mL solution in water
CAS Number 28718-90-3
Molecular Formula $C_{16}H_{17}Cl_2N_5$
Molecular Weight 350.25
Appearance yellow solution
Lot Number AB2040

	LOT DATA	SPECIFICATION
Absorption ¹ Maximum	347 nm	348 ± 3 nm
Fluorescence ¹ Emission Maximum	455 nm	457 ± 4 nm
HPLC Purity	96 % at 254 nm	≥ 95 % at 254 nm
NMR Result	meets specification	consistent with structure

1. Solvent: Methanol.